

KASON 481 Series Inside Release (ISR) Instructions



TEST TO ASSURE LATCH RELEASE FUNCTION. TEST TWICE PER WEEK
Function is critical. Entrapment in freezer can cause severe injury or death.

1

LOW TEMPERATURE GREASE

Ø 19.1mm

LATCH not shown

"L"

WARM SIDE

| Kason Model | Kason Latch | "L" |
|-------------|-------------|---------|
| K481A | K56 | 144.5mm |
| K481B | 0052 | 152.4mm |
| | 0053 | 198.4mm |
| K481C | K58 | 134.1mm |
| | K78, K778 | 143.7mm |

2

| Kason Model | "X" |
|-------------|--------|
| K481A | 60.3mm |
| K481B | 95.3mm |
| K481C | 46mm |

WARM SIDE

90°

Ø 19.1mm

"X"

Door

COLD SIDE

CAUTION!
If trimming to size, a full round end is required for proper activation of latch mechanism.

3

NO SEALANT

3X Ø 5.2mm

COLD SIDE

Door + "X"

| Kason Model | Door Max |
|-------------|----------|
| K481A | 152.4mm |
| K481B | |
| K481C | |